

3D microbiotic skin models to understand and modulate skin colonization and infection

ADEBIOTECH - Romainville | June, 2018

Manon Gault

Care
Creations™

■ - BASF

We create chemistry

Microbiote methods applied to dermocosmetic

Studying skin microbiote and modulation by cosmetic active ingredients require the development of several incremental methods to evaluate:


- Several characteristics of commensal strains as growth, biofilm formation, stress impact such as UV or pollution,
- Interaction between skin cells and germs as adhesion, growth, virulence, inflammation.

CO-CULTURES

Materials & Methods:

- 4 ATCC® strains are amplified using selective media combined to aerobic or anaerobic culture condition: *S. aureus* (6538™), *Staphylococcus epidermidis* (12228™), *S. hominis* (27844™) and *Propionibacterium acnes* (11827™).
- *In tubo* co-cultures of up to 4 species and qPCR quantification after 2 days interaction with or without active ingredients at 37°C .
- Evaluation of relative germs growth and modulations by active ingredient (n=3, T-Test (***) p<0.001, (*) p<0.05).

Strains repartition




Selection of active ingredients favoring the growth of skin benefic bacteria to the detriment of more opportunist bacteria.

BIOFILM FORMATION

Materials & Methods:

- Defined bacterial culture associated to suspension of magnetizable particles are placed in the wells of a microplate incubated for 4-6 h for the bacteria to adhere to the substrate and initiate biofilm formation.
- The microplate is placed on a block containing 96 mini-magnets for one minute before reading.
- Results are expressed as a score of 0 to 20, biofilm present to biofilm absent (BioFilm Ring Test[®] - n = 16, Mann-Whitney Rank Sum Test vs relative control, (***) p<0.001).




When a biofilm is formed, the particles are embedded, no spot is visible and BioFilm Index is close to 0. Quantitative biofilm evaluation is easiest compared to agar based method.

POLLUTION IMPACT ON MICROBIOTE

Materials & Methods:

- *S. epidermidis* or *P. acnes* treated or not by 0.5% active ingredient (AI) during PM2.5 exposure:
 - SE 10⁷ ATCC® 14490™ for 24h in TSB media in aerobia or PA 10⁸ ATCC® 6919™ for 72h in Schaedler media in anaerobia.
 - Particulate matters: 1649b urban dust particles containing mainly polycyclic aromatic hydrocarbons, polychlorinated biphenyls and chlorinated pesticides.
- Bacterial growth assessment by spectrophotometry.

Bacterial growth evaluation after PM exposure



Both germs were impacted by pollution.

Significant protection was obtained only in favor of *S. epidermidis* with the AI.


INFLAMMATION

Materials & Methods:

- **HaCat cell** lines grown until confluence and treated or not for 24h with or without active ingredient in dose effect.
- ***S. aureus*** ($1.5 \cdot 10^7$ ATCC® 12600™) applied for 2h.
- **IL-8 quantification & cell viability** (ELISA & MTT).

IL-8 release after *S. aureus* challenge

Results in % / Control




Statistics:
Mean value \pm SD
n = 6
One way ANOVA
dunnett's method
versus treated control
(***) p<0.001

Within 1 day *S. aureus* induced inflammation.



Significant protection was obtained with a specific hydrolyzed yeast extract.

3D EPIDERMAL MODELS



3D reconstructed microbiotic epidermis

- **Epidermal reconstructions:** 19yo primary keratinocytes grown for 14 days differentiation at the air-liquid interface (ALI) on fibroblast underlaying membrane.
- ***S. epidermidis* alone or in combinaison with *P. acnes* or *S. aureus*** (SE 10^5 ATCC® 12228™ & PA 10^7 ATCC® 11828™ or SA 10^4 CIRI SH1000).
- **Topic treatment** with active ingredient (AIT) or PBS for 1 day.
- **Histology & bacterial growth evaluation.**




Compared to *S. epidermidis*, *S. aureus* already damaged skin surface after 1 day exposure. *S. epidermidis* growth is independent of *P. acnes* and is promoted by the active ingredient.

3D SKIN MODELS

3D reconstructed microbiotic skin

- **Microbiotic Mimeskin®:**
 - 61yo fibroblasts grown on a collagen based matrix for 28 days.
 - 61yo keratinocytes grown on dermal equivalents for 7 days before differentiation at the air-liquid interface for 1 or 2 weeks.
- ***S. epidermidis* and *P. acnes* alone or in combination**
 - SE or PA alone at 10^7 or combined at 10^6 SE + 10^7 PA SE ATCC® 14990™ & PA ATCC® 6919™
 - Seeding after 1 or 2 weeks emersion to mimic wound healing or mature skin.
- **Histology and Bacterial growth evaluation.**

Results with seeding at 1 week post emersion



***S. epidermidis* is mainly at the skin surface whereas *P. acnes* is seen deeper within the skin layers due to its anaerobic growth preference. *P. acnes* did not impact *S. epidermidis* whereas *S. epidermidis* inhibited *P. acnes* growth.**


We create chemistry

Care
Creations™